
Sample Report From Horizon™ DXA System




Advanced Body Composition® Report



Legend

- **1.** Visual image of precise location of bone, lean mass, and fat mass
- 2. Plot of % Body Fat compared to age group
- **3.** Precise Fat Mass, Lean Mass, and % Body Fat measurements of each area:
 - Left Arm
 - Right Arm
 - Torso (Trunk)
 - Abdomen (Android)
 - Hips (Gynoid)
 - Left Leg
 - Right Leg
- **4.** Calculation of Total Mass, Fat Mass, and Lean Mass to give overall % Body Fat
- 5. Estimated amount of visceral fat (the type of fat around internal organs associated with medical disorders such as metabolic syndrome, cardiovascular disease, and type 2 diabetes.

Track % Body Fat, Fat Mass and Lean Mass values over time



Legend

- 1. Graph of % Body Fat values over time
- **2.** % Body Fat values over time
- 3. Fat Mass values over time
- 4. Lean Mass values over time
- **5.** Visual comparison over time of changes in bone, lean mass, and fat mass